

Light Scattering (SLS/DLS)

In addition to the neutron scattering experiments described elsewhere a state-of-the-art light scattering laboratory is available. The use of light scattering is necessary if additional information from another Q and t range is required. The existing experimental equipment allows to carry out static as well as dynamic measurements in the same experiment.

Instrument Details

light source:	Coherent Innova 90/4 Ar ⁺ -ion laser, $\lambda = 488\text{nm}, 514.5\text{nm}$ (max. 1.9W TEM00), etalon mode possible Uniphase HeNe laser, $\lambda = 632.8\text{nm}$ (22mW)
polarizers:	Glan-Thompson prism (B. Halle, Berlin)
goniometer:	ALV/SP125 compact goniometer, $Q = 10^{-1}$
detectors:	fiber optical detection system operating with ALV/SO-SIPD detector systems ("pseudo"-cross correlation) ALV high QE APD (autocorrelation, SLS) EMI 9130/100S B03 (autocorrelation, SLS)
correlator:	ALV-5000E fast version, 319 channel (multi-tau)
time range:	$12.5 \times 10^{-9} - 5.47 \times 10^3$ s
Q-range:	$2.3 \times 10^{-4} - 3.8 \times 10^{-3} \text{ \AA}^{-1}$
Intensity	14kHz (toluene, $\theta=90^\circ$)
coherence factor:	0.92
sample volume:	minimum 250 μl
sample environment:	cryostat 260 ...420K rotating sample holder for nonergodic samples (e.g. gels, glasses)


Instrument Description

In static light scattering (SLS) the integrated intensity of the scattered light is measured as a function of the scattering vector Q , i.e. $S(Q)$. A classical field of application is the characterisation of polymers (radius of gyration R_g , weight average of molecular weight M_w , 2nd virial coefficient A_2).

Dynamic light scattering (DLS), also called quasi-elastic light scattering (QUELS) or photon correlation spectroscopy (PCS) works as a spectroscopic method in time domain; the measured quantity is the autocorrelation function of the scattered light intensity $g^2(Q,t) = \langle I(Q,0)I(Q,t) \rangle / \langle I(Q) \rangle^2$. The large time range (about 10 decades) covered by this method as well as its characteristic length scale ($Q^{-1} \sim 100$ nm) make PCS an important method for the investigation of dynamic processes in the melt and in solution (D_{trans} , D_{rot} , R_h , polydispersity). Typical applications range from biological systems (DNA fragments, microorganisms) to the internal dynamics of polymer networks.


Examples of Measurements

a. SLS




Form factor $P(Q)$ of a hard-sphere colloid, Poly(Methylmethacrylat) in cis-Dekalin, with $R=760\text{nm}$, $10^\circ < \Theta = 154^\circ$, $T=293\text{K}$.

b. DLS



The intensity autocorrelation functions $g^2(Q,t)$ shown here from samples with extremely different relaxation times demonstrate the large time window of dynamical light scattering.

- $\tau=2.5\text{s}$; Polystyrene Latex Standard ($d=304\text{nm}$) in $\text{H}_2\text{O}/\text{Glycerin}$ (23:87), 10^{-4} weight-%, $\Theta=30^\circ$, $T=283\text{K}$.
- $\tau_1=68\text{ms}$, $\tau_2=8.6\mu\text{s}$; 18-arm Polyisoprene star polymer ($M_w=152000$ g/mol) in Methylcyclohexan, $c=0.309\text{g/ml}$, $\Theta=90^\circ$, $T=293\text{K}$
- $\tau=3.3\mu\text{s}$; Polystyrene ($M_w=3300$ g/mol) in Benzene, $c=0.1\text{g/ml}$, $\Theta=150^\circ$, $T=293\text{K}$.



Light scattering set-up